WIDEX **MOMENT**[™]

USER INSTRUCTIONS **WIDEX mRIC CHARGER**

WPT1 model

GB
US
DK
DE
FR
ES
IT
JP
SC
PT

WIDEX[®]
HIGH DEFINITION HEARING

GB	USER INSTRUCTIONS	3
US	USER INSTRUCTIONS	21
DK	BRUGSANVISNING	39
DE	BEDIENUNGSANLEITUNG	57
FR	MODE D'EMPLOI	76
ES	INSTRUCCIONES DE USO	95
IT	LIBRETTO D'ISTRUZIONI	113
JP	取扱説明書	131
SC	用户使用手册	150
PT	INSTRUÇÕES DE UTILIZAÇÃO	168
	FCC AND ISED STATEMENTS	188

IMPORTANT SAFETY INFORMATION

Read these pages carefully before you start using your charger.

Keep the charger and other small parts away from children, pets, and mentally disabled people.

Never expose the charger to temperatures above 50° C (122° F), to prolonged sunlight, open flames or ignition sources.

Place the charger upright on a flat, stable surface when charging and with adequate space around the charger. Do not cover or place near flammable substances.

Never try to open or repair the charger yourself. Contact your hearing care professional if you experience defects.

Use only a WIDEX charger with your hearing aids. Use of a charger not recommended by Widex may compromise the system or result in dangerous situations.

Never leave the charger in direct sunlight and do not immerse it in water.

Never insert other objects than your hearing aids in the charging wells.

The power adapter must comply with IEC 60601-1, IEC 62368-1 or equivalent safety standards acceptable in your region. Always use the supplied USB cable with a power adapter with a USB-A type socket.

YOUR CHARGER

Welcome

Thank you for choosing this product. Your WIDEX mRIC charger is intended to detect and charge the WIDEX MOMENT™ MRR2D hearing aid.


Read these instructions carefully before you start using this device.

NOTE

Your charger may not look exactly like the one illustrated in this booklet. We reserve the right to make any changes we consider necessary.

The charger at a glance

The charger has a charging well for each hearing aid, and two light indicators that show the charging status of the hearing aids.

1. Micro USB port
2. Charging wells
3. Light indicators
4. Hinged cover

Operating, storage and transport conditions

Your charger is designed for indoor use. It is sensitive to extreme conditions such as high temperatures, and it should not be operated, stored or transported in direct sunlight.

Operate the charger within the following conditions:

	Minimum	Maximum
Temperature	0°C (32°F)	30°C (86°F)
Humidity	10% rH	95% rH

Store and transport the charger within the following conditions:


	Minimum	Maximum
Temperature	-20°C (-4°F)	50°C (122°F)
Humidity	10% rH	95% rH
Atmospheric pressure	750 mbar	1060 mbar

You can find technical data sheets and additional information on your charger on <https://global.widex.com>.

OPERATING THE CHARGER

Getting ready for use


The cover

Before use, grab the hinged cover on both sides and lift it upwards.


Then pull it over towards the side with the USB port.


When you are not using the charger, close the hinged cover to protect the charging wells. Simply pull the hinged cover upwards to cover the top before pressing it down.

Power supply


You can connect your charger to different power supplies with 5V DC, 500 mA output:

Use only the included USB cable.


Connect the cable to the charger:

1. Micro USB port

Power on


When the charger is connected to the power supply, the light indicators **flash white 5 times** and then turn off:

The charger is now in standby mode and ready to charge your hearing aid.

Charging the hearing aid

Gently place your hearing aid in the charging well. The hearing aid turns off automatically.

When connected to power, the charger can charge two hearing aids simultaneously; however, each hearing aid can be charged individually, depending on the state of the hearing aid.


This means that if one hearing aid is fully charged and placed in the charger, it will not charge.

However, if the other hearing aid is *not* charged, it will charge.

Charging your hearing aid fully will take approximately 4 hours. If you charge your hearing aid for half an hour, you will have approximately 4 hours of use.


Charging indicators

While the charger is detecting the hearing aid, the light indicator turns **solid yellow**:


Once the charger has detected the hearing aid, the charging will begin.


The light indicator will **flash green** to indicate that the hearing aid is low on battery and less than **40% charged**:

The light indicator will **glow green** to indicate that the hearing aid is more than **40% charged**:

The light indicator will **turn solid green** to indicate that the hearing aid is **100% charged**:

After charging

Avoid removing the hearing aids from the charging wells by pulling the earwire, as this could damage the connection between the earwire and the hearing aid.

Instead, put your fingers round the hearing aid when removing it from the charger.

When you remove the hearing aid from the charger, the hearing aid turns on after a few seconds.

Remember to turn it off if you are not going to use it right away.

Immediately after charging, the hearing aid may temporarily feel warm.

If you are not comfortable with this, wait a few minutes, for the hearing aid to reach room temperature, before putting the hearing aid on.

Overview of light indicators

Light indicators	Status
Flashing white 	The charger is connected to power
Solid yellow 	The charger is detecting the hearing aid
Flashing green 	The hearing aid is charging Low battery, less than 40%
Glowing green 	The hearing aid is charging Battery more than 40%
Solid green 	The hearing aid is fully charged Battery 100% charged

Light indicators**Status**

One indicator flashing red

An error occurred during charging and the charging has stopped
See Troubleshooting

Both indicators flashing red

An error occurred during charging and the charging has stopped
See Troubleshooting

Flashing yellow

The temperature of the hearing aid has become too high during charging and charging is paused
See Troubleshooting

CLEANING

The charger does not generally require cleaning. If necessary, carefully wipe off the charger with a soft dry cloth.

If you notice any dust in the charging well, gently blow it clean. If this is not enough, gently clean the charging well with a cotton bud.

To remove earwax from the charger, carefully wipe with a soft cloth with a bit of mild household detergent.

Never use sharp/pointy objects to clean the charger.

TROUBLESHOOTING

These pages contain advice on what to do if your charger does not work as expected.

Light indicators	Possible cause	Solution
No lights	The charger is not connected	Make sure that the USB cable is connected to a power adapter and then plugged into a wall socket
	The hearing aid is not placed correctly in the charger	Place the hearing aid correctly in the charger
One or both indicators flashing yellow 	The temperature of your surroundings or of the hearing aid exceeds a critical limit and the charging is paused	The charging will continue when the temperature of your surroundings or the hearing aid has dropped. If this happens continually, find a cooler place to charge your hearing aid
One indicator flashing red 	The charger cannot recognise the hearing aid	Remove the incompatible hearing aid

Light indicators	Possible cause	Solution
	The hearing aid has a critical error	Remove the hearing aid, disconnect from the power for 10 seconds and then try again If the charger is still flashing red, contact your hearing care professional
Both indicators flashing red 	The charger has a critical error	Remove the hearing aid, disconnect from the power for 10 seconds and then try again If the charger is still flashing red, contact your hearing care professional

REGULATORY INFORMATION

Intended use

The WIDEX mRIC charger WPT1 model is to be used for charging the WIDEX MOMENT™ MRR2D hearing aid.

EU directives

Directive 2014/53/EU

Hereby, WSAUD A/S declares that the WPT1 is in compliance with the essential requirements and other relevant provisions of Directive 2014/53/EU.

The WPT1 contains a radio transmitter operating at: 6.78 MHz, 13.4 dB μ A/m @ 3 m.

A copy of the Declaration of Conformity according to 2014/53/EU can be found at

<https://global.widex.com/doc>

N26346

Information regarding disposal

Do not dispose of hearing aids, hearing aid accessories and chargers with ordinary household waste.

Hearing aids, hearing aid accessories and chargers should be disposed of at sites intended for electrical waste and electronic equipment, or given to your hearing care professional for safe disposal.

For returning or shipping the devices, please contact your hearing care professional.

SYMBOLS

Symbols commonly used by WSAUD A/S in medical device labelling (labels/IFU/etc.)

Symbol Title/Description

Manufacturer

The product is produced by the manufacturer whose name and address are stated next to the symbol. If appropriate, the date of manufacture may also be stated.

Catalogue number

The product's catalogue (item) number.

Consult instructions for use

The user instructions contain important cautionary information (warnings/precautions) and must be read before using the product.

Warning

Text marked with a warning symbol must be read before using the product.

WEEE mark

“Not for general waste”. When a product is to be discarded, it must be sent to a designated collection point for recycling and recovering to prevent the risk of harm to the environment or human health as a result of the presence of hazardous substances.

Symbol Title/Description

CE mark

The product is in conformity with the requirements set out in European CE marking directives.

RCM mark

The product complies with electrical safety, EMC and radio spectrum regulatory requirements for products supplied to the Australian or New Zealand markets.

FCC logo

The product is compliant with relevant USA 47 CFR Rules.
FCC – Federal Communications Commission.

IMPORTANT SAFETY INFORMATION

Read these pages carefully before you start using your charger.

Keep the charger and other small parts away from children, pets, and mentally disabled people.

Never expose the charger to temperatures above 50° C (122° F), to prolonged sunlight, open flames or ignition sources.

Place the charger upright on a flat, stable surface when charging and with adequate space around the charger. Do not cover or place near flammable substances.

Never try to open or repair the charger yourself. Contact your hearing care professional if you experience defects.

Use only a WIDEX charger with your hearing aids. Use of a charger not recommended by Widex may compromise the system or result in dangerous situations.

Never leave the charger in direct sunlight and do not immerse it in water.

Never insert other objects than your hearing aids in the charging wells.

The power adapter must comply with IEC 60601-1, IEC 62368-1 or equivalent safety standards acceptable in your region. Always use the supplied USB cable with a power adapter with a USB-A type socket.

YOUR CHARGER

Welcome

Thank you for choosing this product. Your WIDEX mRIC charger is intended to detect and charge the WIDEX MOMENT™ MRR2D hearing aid.


Read these instructions carefully before you start using this device.

NOTE

Your charger may not look exactly like the one illustrated in this booklet. We reserve the right to make any changes we consider necessary.

The charger at a glance

The charger has a charging well for each hearing aid, and two light indicators that show the charging status of the hearing aids.

1. Micro USB port
2. Charging wells
3. Light indicators
4. Hinged cover

Operating, storage and transport conditions

Your charger is designed for indoor use. It is sensitive to extreme conditions such as high temperatures, and it should not be operated, stored or transported in direct sunlight.

Operate the charger within the following conditions:

	Minimum	Maximum
Temperature	0°C (32°F)	30°C (86°F)
Humidity	10% rH	95% rH

Store and transport the charger within the following conditions:


	Minimum	Maximum
Temperature	-20°C (-4°F)	50°C (122°F)
Humidity	10% rH	95% rH
Atmospheric pressure	750 mbar	1060 mbar

You can find technical data sheets and additional information on your charger on <https://global.widex.com>.


OPERATING THE CHARGER

Getting ready for use


The cover

Before use, grab the hinged cover on both sides and lift it upwards.


Then pull it over towards the side with the USB port.


When you are not using the charger, close the hinged cover to protect the charging wells. Simply pull the hinged cover upwards to cover the top before pressing it down.

Power supply


You can connect your charger to different power supplies with 5V DC, 500 mA output:

Use only the included USB cable.


Connect the cable to the charger:

1. Micro USB port

Power on


When the charger is connected to the power supply, the light indicators **flash white 5 times** and then turn off:

The charger is now in standby mode and ready to charge your hearing aid.

Charging the hearing aid

Gently place your hearing aid in the charging well. The hearing aid turns off automatically.

When connected to power, the charger can charge two hearing aids simultaneously; however, each hearing aid can be charged individually, depending on the state of the hearing aid.


This means that if one hearing aid is fully charged and placed in the charger, it will not charge.

However, if the other hearing aid is *not* charged, it will charge.

Charging your hearing aid fully will take approximately 4 hours. If you charge your hearing aid for half an hour, you will have approximately 4 hours of use.


Charging indicators

While the charger is detecting the hearing aid, the light indicator turns **solid yellow**:


Once the charger has detected the hearing aid, the charging will begin.


The light indicator will **flash green** to indicate that the hearing aid is low on battery and less than **40% charged**:

The light indicator will **glow green** to indicate that the hearing aid is more than **40% charged**:

The light indicator will **turn solid green** to indicate that the hearing aid is **100% charged**:

After charging

Avoid removing the hearing aids from the charging wells by pulling the earwire, as this could damage the connection between the earwire and the hearing aid.

Instead, put your fingers round the hearing aid when removing it from the charger.

When you remove the hearing aid from the charger, the hearing aid turns on after a few seconds.

Remember to turn it off if you are not going to use it right away.

Immediately after charging, the hearing aid may temporarily feel warm.

If you are not comfortable with this, wait a few minutes, for the hearing aid to reach room temperature, before putting the hearing aid on.

Overview of light indicators

Light indicators

Status

Flashing white

The charger is connected to power

Solid yellow

The charger is detecting the hearing aid

Flashing green

The hearing aid is charging
Low battery, less than 40%

Glowing green

The hearing aid is charging
Battery more than 40%

Solid green

The hearing aid is fully charged
Battery 100% charged

Light indicators**Status**

One indicator flashing red

An error occurred during charging and the charging has stopped
See Troubleshooting

Both indicators flashing red

An error occurred during charging and the charging has stopped
See Troubleshooting

Flashing yellow

The temperature of the hearing aid has become too high during charging and charging is paused
See Troubleshooting

CLEANING

The charger does not generally require cleaning. If necessary, carefully wipe off the charger with a soft dry cloth.

If you notice any dust in the charging well, gently blow it clean. If this is not enough, gently clean the charging well with a cotton bud.

To remove earwax from the charger, carefully wipe with a soft cloth with a bit of mild household detergent.

Never use sharp/pointy objects to clean the charger.

TROUBLESHOOTING

These pages contain advice on what to do if your charger does not work as expected.

Light indicators	Possible cause	Solution
No lights	The charger is not connected	Make sure that the USB cable is connected to a power adapter and then plugged into a wall socket
	The hearing aid is not placed correctly in the charger	Place the hearing aid correctly in the charger
One or both indicators flashing yellow 	The temperature of your surroundings or of the hearing aid exceeds a critical limit and the charging is paused	The charging will continue when the temperature of your surroundings or the hearing aid has dropped. If this happens continually, find a cooler place to charge your hearing aid
One indicator flashing red 	The charger cannot recognise the hearing aid	Remove the incompatible hearing aid

Light indicators	Possible cause	Solution
	The hearing aid has a critical error	Remove the hearing aid, disconnect from the power for 10 seconds and then try again If the charger is still flashing red, contact your hearing care professional
Both indicators flashing red 	The charger has a critical error	Remove the hearing aid, disconnect from the power for 10 seconds and then try again If the charger is still flashing red, contact your hearing care professional

REGULATORY INFORMATION

Intended use

The WIDEX mRIC charger WPT1 model is to be used for charging the WIDEX MOMENT™ MRR2D hearing aid.

EU directives

Directive 2014/53/EU

Hereby, WSAUD A/S declares that the WPT1 is in compliance with the essential requirements and other relevant provisions of Directive 2014/53/EU.

The WPT1 contains a radio transmitter operating at: 6.78 MHz, 13.4 dB μ A/m @ 3 m.

A copy of the Declaration of Conformity according to 2014/53/EU can be found at

<https://global.widex.com/doc>

N26346

Information regarding disposal

Do not dispose of hearing aids, hearing aid accessories and chargers with ordinary household waste.

Hearing aids, hearing aid accessories and chargers should be disposed of at sites intended for electrical waste and electronic equipment, or given to your hearing care professional for safe disposal.

For returning or shipping the devices, please contact your hearing care professional.

Wireless radio transmission

Electromagnetic emissions/immunity standards

The product complies with the following EMC/EMI standards:

Standard	Test type	Note
47 CFR Part 15, Subpart B	RF emissions	US FCC requirements for unintentional radiators
47 CFR Part 18	RF emissions	US FCC requirements for Industrial, Scientific and Medical Equipment
IEC 60601-1-2	General EMC including RF emissions, RF immunity, ESD	Medical Electrical Equipment, Electro-magnetic compatibility

SYMBOLS

Symbols commonly used by WSAUD A/S in medical device labelling (labels/IFU/etc.)

Symbol Title/Description

Manufacturer

The product is produced by the manufacturer whose name and address are stated next to the symbol. If appropriate, the date of manufacture may also be stated.

Catalogue number

The product's catalogue (item) number.

Consult instructions for use

The user instructions contain important cautionary information (warnings/precautions) and must be read before using the product.

Warning

Text marked with a warning symbol must be read before using the product.

WEEE mark

“Not for general waste”. When a product is to be discarded, it must be sent to a designated collection point for recycling and recovering to prevent the risk of harm to the environment or human health as a result of the presence of hazardous substances.

Symbol Title/Description

CE mark

The product is in conformity with the requirements set out in European CE marking directives.

RCM mark

The product complies with electrical safety, EMC and radio spectrum regulatory requirements for products supplied to the Australian or New Zealand markets.

FCC logo

The product is compliant with relevant USA 47 CFR Rules.
FCC – Federal Communications Commission.

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FCC AND ISED STATEMENTS

Federal Communications Commission Statement

This device complies with part 15/part 18 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

NOTE:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.

- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

NOTE:

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

Changes or modifications to the equipment not expressly approved by WSAUD A/S could void the user's authority to operate the equipment.

ISED Statement / Déclaration d'ISED

CAN RSS-216/CNR-216

CAN ICES-1/NMB-1

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

- (1) This device may not cause interference.
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

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ISED RADIATION EXPOSURE STATEMENT:

This equipment complies with ISED RF radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

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